**GB** Attention Operating Surgeon

**IMPORTANT Please read carefully before implanting this product**

#### DESCRIPTION

The Unity Knee™ is a fixed bearing total knee replacement system that consists of a femoral component, tibial insert, tibial tray and all-polyethylene patellar component for use in primary total knee arthroplasty. The Unity Knee™ is provided in two variants, cruciate retaining (CR) and posterior stabilized (PS). The Unity Knee™ CR is intended for use in instances where the posterior cruciate ligament (PCL) is functional. The Unity Knee™ PS variant is utilized when a total knee replacement is indicated, and the posterior cruciate ligament is non-functioning or absent, resulting in joint instability. The Unity Knee™ patellar component is optional and is used in situations where replacement of the articular surface of the patella is required. The system also provides augmant components including femoral augments tibial augments stem extensions, and offset connections.

The Unity Knee™ is intended for use in total knee arthroplasty in skeletally mature patients, to provide increased mobility and reduce pain by replacing the damaged knee joint articulation where there is evidence of sufficient sound bone to

seat and support the components.

### INDICATIONS

- General total knee arthroplasty indications include:
  - Painful, disabling joint disease of the knee resulting from: degenerative arthritis, rheumatoid arthritis or post-traumatic arthritis
  - Post-traumatic loss of knee joint configuration and function
  - Moderate varus, valgus, or flexion deformity in which the ligamentous structures can be returned to adequate function
  - Revision of previous unsuccessful knee replacement or other procedure, where soft tissue stability is adequate
  - Fracture of the distal femur and/or proximal tibia that cannot be stabilized by standard fracture management techniques
  - The posterior stabilized variant is also indicated for PCL instability requiring implant bearing surface geometries with increased anterior-posterior constraint and absent or non-functioning posterior cruciate ligament

The Unity Knee™ is intended for cemented use, single use only.

#### CONTRAINDICATIONS

- Severe instability secondary to advanced loss of osteochondral structure or the absence of collateral ligament integrity
- Infection/ distant foci of infections

A patient with a knee replacement. The knee is the large, rounded area on the right side of the image.

- Osteomyelitis, Osteoporosis, Osteomalacia
- Marked bone loss or bone resorption
- Metabolic disorders which may impair bone formation
- Vascular insufficiency
- Muscular atrophy or neuromuscular disease
- Allergy to implant material
- Severe deformity

#### WARNINGS & PRECAUTIONS

The following conditions tend to adversely affect knee replacement implants: excessive patient weight, high levels of patient activity, likelihood of falls, poor bone stock, metabolic disorders, disabilities of other joints.

Failure of the implant may result from: selection of improper components, mal-alignment or mal-positioning of the components, inadequate fixation of components, failure to remove surgical debris prior to closure. Clean gloves should be worn when handling implants. Avoid scratching or denting implant surfaces. Do not modify implants. Do not reuse implants or disassemble and reassemble implant components, an implant or implant component may appear undamaged but previous stress may have created imperfections that would reduce the service life of the implant. Refer to the surgical technique regarding the use of augment components.

Surgeons should be thoroughly familiar with the Corin Unity Knee™ surgical technique, the

implants, and the instruments prior to performing the surgery. Corin designs specialised instruments for their joint replacement systems to aid in the accurate implantation of the prostheses. The use of instruments or implant components from other systems can result in inaccurate fit, sizing, and device failure. Surgical instruments are subject to wear with normal usage. Instruments which have experienced extensive use or excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose. Components labelled for "Cemented Use Only" are to be implanted only with bone cement. The Corin Unity Knee™ has not been evaluated for safety and compatibility in the MR environment and has not been tested for heating or migration in the MR environment.

#### ADVERSE REACTIONS

Complications can potentially occur following any joint replacement surgery. The following are the most frequent adverse events after knee arthroplasty: change in position of the components, loosening, tibial subsidence, bending, cracking, fracture, deformation or wear of one or more of the components, infection, pain, dislocation, subluxation, flexion contracture, decreased range of motion, looseness or wear of components; fractures of the femur or tibia.

A patient with a knee replacement. The knee is the large, rounded area on the right side of the image.

#### STERILITY

The Unity Knee™ implants are supplied sterile. Check the integrity of the packaging carefully. Do not use if the packaging is open or damaged. Do not resterilise. For single use only. The product is not labelled 'pyrogen free'.

#### INSTRUMENTS

Instruments used to implant the device are supplied non-sterile and must be sterilised prior to use. Instruments should be cleaned according to the Corin cleaning instructions supplied with the instrument set. After cleaning and prior to sterilisation, the reusable device(s) should be double-wrapped or packaged in CSR sterilisation wraps or pouches. Wrapping should be performed using the appropriate wrapping method (e.g. AAMI CSR wrapping technique). The following sterilisation method has been validated, based on AAMI/ANSI/ISO guidelines and recommendations:

- Method: Moist-Heat Sterilisation
- Cycle: Pre-Vacuum (Pre-Vac)
- Temperature: 270 °F (132 °C)
- Exposure Time: 4 minutes
- Pressure: 2-15 PSIA
- Dry-Time: 30 minutes (minimum, in chamber)
- Cool-Time: 60 minutes (minimum, at room temperature)

After sterilising, the reusable device(s) should be left on the steriliser cart, untouched, for a

A patient with a knee replacement. The knee is the large, rounded area on the right side of the image.

minimum of one hour at room temperature, or until adequately cooled for safe handling. Other sterilisation methods and cycles may also be suitable, however, individuals or hospitals are advised to validate whichever method they deem appropriate at their institution. Ethylene Oxide (EtO) sterilisation and cold sterilisation techniques are not recommended.

The information contained in this package insert was current on the date it was printed but the package insert may have been revised after that date. To obtain a current package insert, please contact Corin.

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