GB Attention Operating Surgeon

3 Biolox Option modular heads are not approved for use in the USA.

IMPORTANT Please read carefully before implanting this product

DESCRIPTION

The Trinity Acetabular System is a modular acetabular cup system consisting of a press fit, titanium alloy shell for use with either ceramic¹ or polyethylene liners and a dedicated range of ceramic (i.e. Biolox Delta) and cast cobalt chrome alloy modular 12/14 taper femoral heads providing ceramic on ceramic¹, ceramic on polyethylene and metal on polyethylene articulations for use in total hip replacement procedures using any Corin femoral stem² with a 12/14 taper connection. The acetabular shell is coated with a rough titanium plasma spray with an additional top layer of advanced biomimetic coating. The acetabular shell is provided with screw holes permitting the use of dedicated titanium bone screws to provide additional fixation if required. Titanium occluders are provided to occlude unused screw holes and an apical introducer hole, where appropriate. The Trinity Acetabular System is intended for use in total hip arthroplasty in skeletally mature patients, to provide increased mobility and reduce pain by replacing the damaged hip joint articulation where there is evidence of sufficient sound bone to seat and support the components.

1 Ceramic liners and ceramic on ceramic articulations are not approved for use in the USA.

2 Ceramic femoral heads are indicated only with Corin titanium stems in the USA.

INDICATIONS

The Trinity Acetabular System is indicated for use

in non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis, rheumatoid arthritis, correction of functional deformity, developmental dysplasia of the hip (DDH) or congenital dysplasia of the hip (CDH). The Trinity Acetabular System is intended for cementless, single use only.

CONTRAINDICATIONS

- Infection
- Osteomyelitis
- Sepsis
- Osteomalacia
- Distant foci of infections
- Osteoporosis
- Marked bone loss or bone resorption
- Metabolic disorders which may impair bone formation
- Vascular insufficiency
- Muscular atrophy or neuromuscular disease
- Allergy to implant material
- Severe deformity

WARNINGS AND PRECAUTIONS

Failure of the implant may result from: selection of improper components, mal-alignment or mal-positioning of the components, inadequate fixation of components, failure to remove surgical debris prior to closure. Clean gloves should be worn when handling implants. Avoid scratching or denting implant surfaces. Do not modify implants. Do not reuse implants, an implant may appear undamaged but previous stress may have created imperfections that would reduce the service life of the implant.

Patients should be warned that ceramic implants can break or become damaged as a result of strenuous activity or trauma, etc. Only the

compatible Corin 12/14 femoral stems components identified in the device description section should be used with the Corin ceramic femoral heads. The use of other femoral stem components may result in fracture or damage to the ceramic femoral head and lead to early failure of the system. Carefully examine each component and its packaging for any signs of damage that may have occurred during shipping or handling (e.g. if dropped on the floor or if scratched by an instrument). Do not implant components if the packaging is damaged or if the implant shows signs of damage. The use of damaged components may lead to premature failure of the device. Be careful when handling ceramic components during assembly; due to the brittle nature of the material, ceramic components are particularly susceptible to premature failure when scratched, cracked or otherwise damaged. Ensure that surfaces of the implant that lock together are clean and dry during implantation in order to ensure proper seating and assembly. Once the femoral stem taper has been assembled to a ceramic head, it should not be reassembled to another ceramic head. If the ceramic head is chipped, cracked, or damaged during head/ stem assembly, the physician should replace both the ceramic head and the femoral stem. In the case of revision due to a fractured ceramic component, when keeping a well fixed stem, the modular head must be replaced with a Biolox Option³ modular head where available. If a Biolox Option modular head is not available, in addition to revising the modular head, the femoral stem must also be revised to provide a suitable taper to which a new ceramic femoral modular head can be attached. In the event of revision not due to fracture of a ceramic component, when keeping a well fixed stem, the modular head must be replaced with a metal modular head or a Biolox Option modular head only.

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